



# CANNON<sup>®</sup> APPLICATION NOTE #0609291

---

**ISSUE: CFR Compliance**

**DATE: 9/29/2006**

---

**Applicable instrument(s): Various**

**ASTM ref. N/A**

---

This Application Note provides information to customers interested in ensuring CANNON software compatibility with their CFR 21 (Part 11) compliance goals.

CANNON Instrument Company is committed to providing software solutions for customers seeking GxP, CFR 21 remediation and full compliance. The following statement describes CANNON software with respect to CFR compliance issues:

- ✓ VISCPRO<sup>®</sup> controlling software for CANNON instruments offers three security levels with login accountability and administrator-configurable accounts for individual users.
- ✓ All data is contained in a Microsoft<sup>®</sup> Access<sup>®</sup>-type database that can be controlled, archived and monitored through the customer quality system. This can even include running VISCPRO from a directory in which only administrators have write privileges.
- ✓ VISCPRO runs independently of the Microsoft Access software application, making it easy for administrators to control rights to the database by limiting the number of applications installed or accessible from the host computer.
- ✓ VISCPRO does not allow modification of measured data by the VISCPRO application.
- ✓ Time-stamped paper copy may be generated for any data set, permitting hard copy archiving and verification of database integrity.
- ✓ A separate Database Manager program provides convenient archiving/retention of data sets.

Specific questions regarding CANNON software and related CFR compliance issues may be addressed to:

**Quality Control Director  
CANNON Instrument Company  
2139 High Tech Road  
State College, PA 16803**