



2139 High Tech Road • State College • PA • 16803 • USA  
814 353 8000 • 800 676 6232 • FAX 814 353 8007

E-MAIL: [cannon@cannoninstrument.com](mailto:cannon@cannoninstrument.com) • WEB: [www.cannoninstrument.com](http://www.cannoninstrument.com)

Date:	December 8, 2006	Rev:	A
CONTACT:	Frank Saylor	Phone:	814-353-8000
	Vice President of Sales and Marketing	Fax:	814-353-8007
		WEB:	<a href="http://www.cannoninstrument.com">www.cannoninstrument.com</a>
		E-MAIL:	<a href="mailto:cannon@cannoninstrument.com">cannon@cannoninstrument.com</a>

RE: **Energy Saver**

## **CAV-2100/2200 Energy Saver increases efficiency**

Available as of 1/1/07 the VISCPRO II for Windows software application will provide a unique “Energy Saver” option that permits the user to configure the CAV-2100/2200 automatic viscometer for automatic shutdown and restart of selected primary systems, including bath heaters, ballast, and instrument display.

A seven-day programmable interface allows for setting ON/OFF times for the instrument. Smart technology ensures that the instrument engages required systems at the optimum time to ensure that the CAV-2100/2200 viscometer is ready to test samples at the desired test temperature at the desired ON time.

The CAV onboard microprocessor will ensure that regardless of computer on/off status, the instrument will maintain the programmed schedule.

###