



2139 High Tech Road • State College • PA • 16803 • USA
814 353 8000 • 800 676 6232 • FAX 814 353 8007

E-MAIL: cannon@cannoninstrument.com • WEB: www.cannoninstrument.com

Date:	July 15, 2008	Rev:	A
CONTACT:	Frank Saylor	Phone:	814-353-8000
	Vice President of Sales and Marketing	Fax:	814-353-8007
		WEB:	www.cannoninstrument.com
		E-MAIL:	cannon@cannoninstrument.com

RE: **PRESS RELEASE—For Immediate Release**

CANNON introduces miniQV™-X quick automatic viscometer

A new approach to in-service oil viscosity measurement



CANNON Instrument Company is excited to announce the release of a quick fully automatic viscometer, the miniQV-X. The new instrument stands just over 20 inches tall and has a removable 25 position rotary auto sampler. Occupying roughly the same footprint as a tabletop rotational viscometer, the miniQV-X automates the time-consuming sample measurement and viscometer tube wash/dry procedures associated with the ASTM D445 method.

Operation of the miniQV-X is simple and quick with typical sample-to-sample cycle times of less than three minutes. This instrument has been developed for the used oil analysis market and offers customers an alternative to slower but fully compliant ASTM D445 instruments manufactured by CANNON.