



2139 High Tech Road • State College • PA • 16803 • USA
814 353 8000 • 800 676 6232 • FAX 814 353 8007

E-MAIL: cannon@cannoninstrument.com • WEB: www.cannoninstrument.com

Date:	June 17, 2008	Rev:	A
CONTACT:	Frank Saylor	Phone:	814-353-8000
	Vice President of Sales and Marketing	Fax:	814-353-8007
		WEB:	www.cannoninstrument.com
		E-MAIL:	cannon@cannoninstrument.com

RE: **PRESS RELEASE—For Immediate Release**

CANNON updates miniAV Series automatic viscometer



CANNON Instrument Company has announced the release of a new generation fully automatic viscometer, the miniAV-X. The product is the newest in the miniAV Series. The engineering update to the original miniAV instrument provides sample-handling capability for multiple-sample testing. The carousel accepts 10 vials and only requires five mL of sample per test. As with the miniAV, the miniAV-X™, automates sample measurements and tube cleaning, the most labor-intensive parts of the kinematic viscosity measurement process.