

Section 1. Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product Name FPRM9D

1.2. Relevant identified uses of the substance or mixture and uses advised against

Recommended Use Flash point testing equipment performance verification reference material

Uses advised against No information available

1.3. Details of the supplier of the safety data sheet

Supplier

Cannon Instrument Company
2139 High Tech Rd.
State College, PA 16803-1733
TEL: (814) 353-8000; (800) 676-6232

For further information, please contact

E-mail Address No information available.

1.4. Emergency telephone number

Emergency Telephone Number (800) 255-3924 Domestic CHEM-TEL Inc.
+1 (813) 248-0585 Overseas CHEM-TEL Inc. (Please Call Collect)

Europe	112
--------	-----

Section 2. Hazards identification

2.1. - Classification of the substance or mixture

REGULATION (EC) No 1272/2008

Chronic Aquatic Toxicity	Category 4
--------------------------	------------

Physical Hazards

None

2.2. Label Elements

Hazard Statements

H413 - May cause long lasting harmful effects to aquatic life

Precautionary Statements

P273 - Avoid release to the environment
P501 - Dispose of contents/ container to an approved waste disposal plant

2.3. Other information

May be irritating to eyes, respiratory system and skin.

Section 3. Composition/information on ingredients

3.1. Substances

3.2. Mixtures

Chemical Name	EC-No	CAS-No	Weight %	EU - GHS Substance	REACH No.
---------------	-------	--------	----------	--------------------	-----------

				Classification	
1-Dodecene, trimer, hydrogenated	417-070-7	151006-62-1	50-80	Aquatic Chronic 4 (H413)	No data available
1-Dodecene, homopolymer, hydrogenated	-	151006-63-2	30-50		No data available

For the full text of the H-Statements mentioned in this Section, see Section 16

Section 4. First aid measures

4.1. Description of first-aid measures

Eye Contact	Rinse thoroughly with plenty of water, also under the eyelids. Get medical attention if irritation persists.
Skin Contact	Wash off with warm water and soap. Get medical attention if symptoms occur.
Ingestion	Do NOT induce vomiting. Drink plenty of water. Never give anything by mouth to an unconscious person. Get medical attention.
Inhalation	IF INHALED: If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing. Get medical attention.

4.2. Most important symptoms and effects, both acute and delayed

Most Important Symptoms/Effects Irritation.

4.3. Indication of immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

Section 5. Fire-fighting measures

5.1. Extinguishing media

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Extinguishing media which must not be used for safety reasons

No information available.

5.2. Special hazards arising from the substance or mixture

Special exposure hazards arising from the substance or preparation itself, combustion products, resulting gases

None in particular.

5.3. Advice for firefighters

Special protective equipment for fire-fighters

As in any fire, wear self-contained breathing apparatus and full protective gear.

Section 6. Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing. Refer to Section 8 for personal protective equipment.

6.2. Environmental precautions

Prevent further leakage or spillage if safe to do so.

6.3. Methods and materials for containment and cleaning up

Prevent further leakage or spillage if safe to do so.

Dam up. Contain and collect spillage with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and place in container for disposal according to local/national regulations (see Section 13). Clean contaminated surface thoroughly.

6.4. Reference to other sections

See Section 12 for additional information.

Section 7. Handling and storage

7.1. Precautions for Safe Handling

Handling

Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin, eyes and clothing. Remove and wash contaminated clothing before re-use. Wear personal protective equipment.

Hygiene Measures

Handle in accordance with good industrial hygiene and safety practice.

7.2. Conditions for safe storage, including any incompatibilities

Keep containers tightly closed in a dry, cool and well-ventilated place.

7.3. Specific end use(s)

Exposure Scenario

No information available.

Other Guidelines

No information available.

Section 8. Exposure controls/personal protection

8.1. Control parameters

Exposure Limits

This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

Biological occupational exposure limits

This product, as supplied, does not contain any hazardous materials with biological limits established by the region specific regulatory bodies

Derived No Effect Level

No information available.

Predicted No Effect Concentration (PNEC)

No information available.

8.2. Exposure controls

Engineering Measures

Ensure adequate ventilation, especially in confined areas.

Personal protective equipment

Personal protection equipment should be chosen according to the CEN standards

Eye Protection

Safety glasses with side-shields.

Skin and Body Protection

Long sleeved clothing.

Hand Protection

Protective gloves.

Respiratory Protection

When workers are facing concentrations above the exposure limit they must use

appropriate certified respirators.

Environmental Exposure Controls No information available.

Section 9. Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical State	Liquid	Appearance	Colorless
Odor	Faint hydrocarbon		

<u>Property</u>	<u>Values</u>	<u>Remarks/ - Method</u>
pH	No data available	None known
Melting Point/Range	No data available	None known
Boiling Point/Boiling Range	>260 °C	None known
Flash Point	> 246 °C	Open cup
Evaporation rate	No data available	None known
Flammability (solid, gas)	No data available	None known
Flammability Limits in Air	No data available	None known

Vapor Pressure	No data available.	None known
Vapor Density	No data available.	None known
Relative Density	No data available 0.84	None known
Water Solubility	Insoluble in water.	None known
Solubility in other solvents	Soluble in solvents.	None known
Partition coefficient: n-octanol/water	No data available	None known
Autoignition Temperature	No data available	None known
Decomposition Temperature	No data available	None known
Viscosity	38 cSt @ 40°C	None known

Explosive Properties	No information available
Oxidizing Properties	No information available

9.2. Other information

VOC Content (%)	No information available
------------------------	--------------------------

Section 10. Stability and reactivity

10.1. Reactivity

Not reactive under normal conditions.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

10.4. Conditions to avoid

Keep away from open flames, hot surfaces and sources of ignition.

10.5. Incompatible materials

Strong oxidizing agents.

10.6. Hazardous decomposition products

Carbon oxides.

Section 11. Toxicological information

11.1. Information on toxicological effects

Acute Toxicity

Product Information

Inhalation

Inhalation of mist may cause irritation to the respiratory system.

Eye Contact

Contact with eyes may cause irritation.

Skin Contact

May cause skin irritation and/or dermatitis. Repeated or prolonged contact may cause

Ingestion	localized dermal effects including contact dermatitis, dry skin, or rash. Ingestion may cause gastrointestinal irritation, nausea, vomiting and diarrhea.
Acute Toxicity	180% of the mixture consists of ingredient(s) of unknown toxicity.

The following values are calculated based on chapter 3.1 of the GHS document:

LD50 Oral	99,999.00 mg/kg
LD50 Dermal	99,999.00 mg/kg
Gas	99,999.00 mg/L
Dust/Mist	99,999.00 mg/L
Vapor	99,999.00 mg/L

Sensitization	No information available.
Mutagenic Effects	No information available.
Carcinogenic Effects	No information available.
Reproductive Toxicity	No information available.
Developmental Toxicity	No information available.
STOT - single exposure	No information available.
STOT - repeated exposure	No information available.
Aspiration Hazard	No information available.

Section 12. Ecological information

12.1. Toxicity

Ecotoxicity Effects

Contains no substances known to be hazardous to the environment at concentrations that would be significant.

12.2. Persistence and degradability

No information available.

12.3. Bioaccumulative potential

No information available.

12.4. Mobility in soil

No information available.

12.5. Results of PBT and vPvB assessment

No information available.

12.6. Other adverse effects

This product does not contain any known or suspected endocrine disruptors.

Section 13. Disposal considerations

13.1. Waste treatment methods

Waste from Residues / Unused Products

Dispose of in accordance with local regulations.

Contaminated Packaging

Empty containers should be taken to an approved waste handling site for recycling or

disposal.

Section 14. Transport information

Note:

The information provided below may not apply to all shipping situations. Consult appropriate Dangerous Goods Regulations for additional requirements and mode-specific, material-specific, or quantity-specific shipping requirements.

IMDG/IMO

14.1. UN-Number	Not regulated.
14.2. Proper Shipping Name	Not regulated.
14.3. Hazard Class	Not regulated.
14.4. Packing Group	Not regulated.
Description	Not applicable.
14.5. Marine Pollutant	None
14.6. Special Provisions	None
14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	No information available.

RID

14.1. UN-Number	Not regulated.
14.2. Proper Shipping Name	Not regulated.
14.3. Hazard Class	Not regulated.
14.4. Packing Group	Not regulated.
Description	Not applicable.
14.5. Environmental hazard	None
14.6. Special Provisions	None

ADR

14.1. UN-Number	Not regulated.
14.2. Proper Shipping Name	Not regulated.
14.3. Hazard Class	Not regulated.
14.4. Packing Group	Not regulated.
Description	Not applicable.
14.5. Environmental hazard	None
14.6. Special Provisions	None

ICAO

14.1. UN-Number	Not regulated.
14.2. Proper shipping name	Not regulated.
14.3. Hazard Class	Not regulated.
14.4. Packing Group	Not regulated.
Description	Not applicable.
14.5. Environmental hazard	None
14.6. Special Provisions	None

IATA

14.1. UN-Number	Not regulated.
14.2. Proper Shipping Name	Not regulated.
14.3. Hazard Class	Not regulated.
14.4. Packing Group	Not regulated.
Description	Not applicable.
14.5. Environmental hazard	None
14.6. Special Provisions	None

Section 15. Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**International Inventories**

TSCA	Complies
EINECS/ELINCS	Complies

DSL/NDSL	Complies
PICCS	Complies
ENCS	Complies
IECSC	Complies
AICS	Complies
KECL	Complies

Legend

TSCA - United States Toxic Substances Control Act Section 8(b) Inventory

EINECS/ELINCS - European Inventory of Existing Commercial Chemical Substances/EU List of Notified Chemical Substances

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

PICCS - Philippines Inventory of Chemicals and Chemical Substances

ENCS - Japan Existing and New Chemical Substances

IECSC - China Inventory of Existing Chemical Substances

AICS - Australian Inventory of Chemical Substances

KECL - Korean Existing and Evaluated Chemical Substances

15.2. Chemical Safety Assessment

No information available

Section 16. Other information**Full text of H-Statements referred to under sections 2 and 3**

H413 - May cause long lasting harmful effects to aquatic life

Key literature references and sources for data

www.ChemADVISOR.com/

Issuing Date	24-Apr-2014
Revision Date	29-Aug-2017
Revision Note	(M)SDS sections updated: 3, 8.

This safety data sheet complies with the requirements of Commission Regulation (EU) 2015/830 of 28 May 2015 amending Regulation (EC) No. 1907/2006

General Disclaimer

The information provided on this SDS is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guide for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered as a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

End of Safety Data Sheet